



Measuring chart

Spinomed® active
Spinomed® active men

measuring chart Spinomed® active

Please measure while patient is wearing a bra!




below bust
horizontally
just below the breast



bust
horizontally over the
widest part of the breast




bust horizontally
from one side of the
breast cup over the
middle to the end of the
breast



bust vertically
from below bust to the
middle



shoulder bust
close fitting from middle
of shoulder to the middle
of bust



waistline
circumference of waist




hip
circumference at the
widest part of the hip

Tips for measuring


- please measure below the bust so that the patient is able to breath easily
- for exact measurement of patients with wide bust size please ease the bra supports
- the arms hang freely to the sides
- please indicate if patient was not wearing a bra
- in any case of doubt, please send photos



frontal height
below bust to thigh
crease



length of the back 1
from C7 to sacrum
following the contours
of the back



length of the back 2
close fitting from C7 to
the widest part of the
hip, then hanging
vertically to the base of
buttocks

Spinomed® active


customer _____
 address _____
 doctor _____
 insurance company _____
 measured by _____ Stamp _____

medi GmbH & Co. KG
 Medicusstraße 1
 95448 Bayreuth
 Germany
 T +49 921 912-0
 F +49 921 912-783
 export@medi.de
 www.medi.de

Please measure while patient is wearing a bra


For Spinomed active without cup only the measurement in the grey are necessary

with without gusset: hooks/eyes velcro/coating first treatment further treatment without backbrace




below bust horizontally
just below the breast

_____ cm



bust horizontally
horizontally over the widest part of the breast

_____ cm




bust horizontally
from one side of the breast cup over the middle to the end of the breast

_____ cm




bust vertically
from below bust to the middle

_____ cm




shoulder bust
close fit from middle of shoulder to middle of bust

_____ cm




waistline
circumference of waist

_____ cm




hip circumference
at the widest part of the hip

_____ cm




frontal height
below bust to thigh crease

_____ cm



length of the back I
from C7 to sacrum following the contours of back

_____ cm



length of the back II
close fitting from C7 to the widest part of the hip, then hanging vertically to the base of buttocks

_____ cm

normal back normal belly normal bottom patient's height _____ cm
 kyphosis sagging belly flat bottom patient's age _____ years
 hollow back pot belly large bottom patient's weight _____ kg
 size of bra _____

Commission (Please write legibly) _____ (Patients data will be kept confidential under the guidelines of data security)

Comments: (e. g. additional disease) _____

date _____ signature/seal _____

Spinomed® active men

customer _____

address _____

doctor _____


insurance company _____

measured by _____

Stamp _____


medi GmbH & Co. KG
 Medicusstraße 1
 95448 Bayreuth
 Germany
 T +49 921 912-0
 F +49 921 912-783
 export@medi.de
 www.medi.de

first treatment further treatment
without backbrace



bust circumference
of the widest part
of the chest

_____ cm




waistline circumference
of waist

_____ cm




hip circumference
of the widest part
of the hip

_____ cm




back I
C7 to tailbone close
fitting

_____ cm




back II
close fitting from C7 to the widest
part of the hip, then hanging
vertically to the base of buttocks

_____ cm




**length from middle of
shoulder to waist**
close fitting
measurement

_____ cm



**length from waist
to pubic bone**
close fitting
measurement

_____ cm



circumference of thigh
please measure the widest
part of thigh

_____ cm

normal back normal belly normal bottom patient's height _____ cm

kyphosis sagging belly flat bottom patient's age _____ years

hollow back pot belly large bottom patient's weight _____ kg

Commission (Please write legibly) _____ (Patients data will be kept confidential under the guidelines of data security)

Comments: (e. g. additional disease) _____

date _____ signature/seal _____

medi. I feel better.